

MERISOFT CHROMIC™ / MERISOFT CHROMIC™ CAS

Chromic Catgut Suture

Properties

MERISOFT CHROMIC™ gut suture is treated with a chromium salt solution to resist the enzymatic action, tensile strength retains until 10 - 14 days and absorption between 90-100 days.

Indications


MERISOFT CHROMIC™ suture is indicated for use in general soft tissue approximation and/or ligation.

Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability



Product Overview

Primary packaging	: Packed in an aluminum foil pouch Non-needled suture packed in a cassette
Secondary packaging	: Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film
Sizes	: USP 6/0 (metric 1) - USP 2 (metric 6)
Colour	: Brown dyed 
Sterilization	: Gamma radiation