

# FILAPRON<sup>TM</sup> VET

## Poliglecaprone 25 Suture

### Properties

- For dyed, 60 - 80% tensile strength retained until 7 days and 30 - 40% retained until 14 days
- For undyed 50 - 60% retained until 7 days and 20 - 30% retained until 14 days.
- Complete mass absorption within 90 to 120 days

### Indications

When short term wound support is desirable


- Skin closure (intracutaneous, subcutaneous)
- Ligatures
- Soft Tissues

### Advantages

- Short-term absorption profile
- High knot pull tensile strength
- Superior handling features



### Product Overview

Presentation/Format	: Dozen suture packed in a Dozen box
Structure	: Monofilament
Colour	: Dyed, Undyed
	
Chemical composition	: Poly(glycolide-co-caprolactone)PGCL
Sizes	: USP 6/0 (metric 0.7) - USP 1 (metric 4)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide